


## MMP2 (3A8) Mouse mAb

ABZ4267


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | MMP2 (3A8) Mouse mAb                       |
| <b>Accession(Primary)</b>   | P08253                                     |
| <b>Host Species</b>         | Mouse                                      |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG1                                       |
| <b>Species Reactivity</b>   | Human, Mouse                               |
| <b>Calculated MW</b>        | Calculated MW: 74 kDa; Observed MW: 74 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IF-1:50-1:200             |

### Target information (P08253)

**Gene Name** MMP2**About Gene Symbol**

Western blot analysis of MMP2 (3A8) in Raw264.7 lysates using MMP2 antibody.



Immunofluorescence analysis of MMP2 (3A8) in Hela using MMP2 antibody, and DAPI (blue).