


## MSRB3 Rabbit mAb

ABZ1998


### Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | MSRB3 Rabbit mAb                           |
| <b>Accession(Primary)</b>   | Q8IXL7                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 21 kDa; Observed MW: 21 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                            |

### Target information (Q8IXL7)

**Gene Name** MSRB3**About Gene Symbol**

Western blot analysis of MSRB3 in 3T3 lysates using MSRB3 antibody.



Western blot analysis of MSRB3 in C6 lysates  
using MSRB3 antibody.