

## MST4 Rabbit mAb


ABZ1786

### Product Overview


|                             |                                                                                          |
|-----------------------------|------------------------------------------------------------------------------------------|
| <b>Name</b>                 | MST4 Rabbit mAb                                                                          |
| <b>Accession(Primary)</b>   | Q9P289                                                                                   |
| <b>Host Species</b>         | Rabbit                                                                                   |
| <b>Clonality</b>            | Monoclonal                                                                               |
| <b>Isotype</b>              | IgG                                                                                      |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                                        |
| <b>Calculated MW</b>        | Calculated MW: 47 kDa; Observed MW: 52 kDa                                               |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:50-1:100;<br>ICC/IF-1:200-1:500; FC-1:100-1:1000; IP-1:20-1:50 |

### Target information (Q9P289)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | STK26 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of MST4 in K562, C6 lysates using MST4 antibody.



Immunocytochemistry analysis of MST4(green) in Hela using MST4 antibody, and DAPI(blue)