

## MUC4 Rabbit mAb


ABZ1891

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | MUC4 Rabbit mAb                              |
| <b>Accession(Primary)</b>   | Q99102                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 542 kDa; Observed MW: 120 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000                             |

### Target information (Q99102)

|                          |                                                                                                                                                                                                                                                                                  |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MUC4                                                                                                                                                                                                                                                                             |
| <b>About Gene Symbol</b> | MUC4 (Mucin 4, Cell Surface Associated) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of MUC4 in PC-3 lysates using Mucin 4 antibody.