


## MYH6/7 Rabbit pAb

ABZ1559


### Product Overview

|                             |                                                         |
|-----------------------------|---------------------------------------------------------|
| <b>Name</b>                 | MYH6/7 Rabbit pAb                                       |
| <b>Accession(Primary)</b>   | P13533                                                  |
| <b>Host Species</b>         | Rabbit                                                  |
| <b>Clonality</b>            | Polyclonal                                              |
| <b>Isotype</b>              | IgG                                                     |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                       |
| <b>Recommended Dilution</b> | IHC-1:50-1:200; ELISA-1:10000-1:20000;<br>IF-1:50-1:200 |


### Target information (P13533)

**Gene Name** MYH6/7**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human heart using MYH6/7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using MYH6/7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human skeletal muscle using MYH6/7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.