

# Mannose Receptor/CD206 Rabbit mAb


ABZ3044

## Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Mannose Receptor/CD206 Rabbit mAb            |
| <b>Accession(Primary)</b>   | P22897                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 166 kDa; Observed MW: 166 kDa |
| <b>Recommended Dilution</b> | WB-1:5000-1:50000                            |

## Target information (P22897)

|                          |                                                                                                                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | MRC1                                                                                                                                                                                                                                                                      |
| <b>About Gene Symbol</b> | MRC1 (Mannose Receptor C-Type 1) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Mannose Receptor in Jurkat, C6, HeLa lysates using Mannose Receptor antibody.