

Myoblast Determination Protein 1 Rabbit pAb

ABZ1669

Product Overview

Name	Myoblast Determination Protein 1 Rabbit pAb
Accession(Primary)	P15172
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat, Chicken
Calculated MW	Calculated MW: 35 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P15172)

Gene Name MYOD1

About Gene Symbol



