

N4BP1 Rabbit mAb

ABZ2029


Product Overview

Name	N4BP1 Rabbit mAb
Accession(Primary)	Q6A037
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 99 kDa; Observed MW: 99 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500

Target information (Q6A037)

Gene Name N4Bp1

About Gene Symbol



Western blot analysis of N4BP1 in K562, HeLa lysates using N4BP1 antibody.

