

NDUFAF1 Rabbit mAb


ABZ1734

Product Overview

Name	NDUFAF1 Rabbit mAb
Accession(Primary)	Q9Y375
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 38 kDa; Observed MW: 35 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:500; FC-1:100-1:200

Target information (Q9Y375)

Gene Name	NDUFAF1
About Gene Symbol	



Western blot analysis of NDUFAF1 in rat Brain, MCF-7, C6, 3T3, HeLa lysates using NDUFAF1 antibody.