

NODAL Rabbit mAb

ABZ1895


Product Overview

Name	NODAL Rabbit mAb
Accession(Primary)	Q96S42
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:1000-1:5000; FC-1:100-1:500

Target information (Q96S42)

Gene Name NODAL

About Gene Symbol



Western blot analysis of Nodal in C6, 3T3, HeLa lysates using Nodal antibody.

