


NOP10 Rabbit mAb

ABZ1822


Product Overview

Name	NOP10 Rabbit mAb
Accession(Primary)	Q9NPE3
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 8 kDa; Observed MW: 8 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:200-1:500

Target information (Q9NPE3)

Gene Name NOP10**About Gene Symbol**

Western blot analysis of NOP10 in HeLa, HL-60 lysates using NOP10 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using NOP10 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.