

NT 3 Human

Neurotrophin-3 Human Recombinant
NTR0062

Product Overview

Name NT 3 Human

Description

Neurotrophin-3 Human Recombinant

Accession (Primary) [P20783](#)

Synonyms

Neurotrophic factor, Nerve growth factor-2, NGF-2, HDNF, NT-3, Neurotrophin-3, Ntf3, Ntf-3, AI316846, AI835689, Nt3.

Introduction

NT3 a member of the neurotrophin family, that controls survival and differentiation of mammalian neurons. This protein is closely related to both nerve growth factor and brain-derived neurotrophic factor. It may be involved in the maintenance of the adult nervous system, and may affect development of neurons in the embryo when it is expressed in human placenta. NTF3-deficient mice generated by gene targeting display severe movement defects of the limbs. The mature peptide of this protein is identical in all mammals examined including human, pig, rat and mouse.

Source

Escherichia Coli.

Physical Appearance

Sterile Filtered White lyophilized (freeze-dried) powder.

Formulation

Lyophilized from 0.02% TFA.

Stability

Lyophilized NGF2 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution NGF-2 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Purity

Greater than 97.0% as determined by (a) Analysis by RP-HPLC. (b) Analysis by SDS-PAGE.

Amino acid sequence

YAEHKSHRGE YSVCDSESLW VTDKSSAIDI RGHQVTVLGE IKTGNPQVKQ YFYETRCKEA RPVKNGCRGI
DDKHWNSQCK TSQTYVRALT SENNKLVGWR WIRIDTSCVC ALSRKIGRT.

Biological Activity

The activity, as determined by the dose-dependent proliferation of BaF3 cells transfected with the TrkB receptor, is typically in the range of 1-10 ng/ml, corresponding to a specific activity of 100,000-1,000,000 units/mg.

Solubility

It is recommended to reconstitute the lyophilized Neurotrophin-3 in sterile 18M²-cm H₂O not less than 100 µg/ml, which can then be further diluted to other aqueous solutions.

Precautions

NT 3 Human is for research use only and not for use in diagnostic or therapeutic procedures.

Target Information: ([P20783](#))