

PEDF Human, His

Pigment Epithelium-Derived Factor Human Recombinant, His Tag

NTR0069

Product Overview

Name	PEDF Human, His
Accession(Primary)	P36955
Immunogen	Anti-human PEDF mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human PEDF protein 20-418 amino acids purified from E. coli.
Description	Pigment Epithelium-Derived Factor Human Recombinant, His Tag

Target information (P36955)

Gene Name

About Gene Symbol