


## NUT Rabbit pAb

ABZ1649

### Product Overview

|                             |                                             |
|-----------------------------|---------------------------------------------|
| <b>Name</b>                 | NUT Rabbit pAb                              |
| <b>Accession(Primary)</b>   | Q86Y26                                      |
| <b>Host Species</b>         | Rabbit                                      |
| <b>Clonality</b>            | Polyclonal                                  |
| <b>Isotype</b>              | IgG                                         |
| <b>Species Reactivity</b>   | Human, Mouse                                |
| <b>Calculated MW</b>        | Calculated MW: 120 kDa; Observed MW: 63 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; ELISA-1:10000-1:20000      |

### Target information (Q86Y26)

**Gene Name** NUTM1**About Gene Symbol**

Western blot analysis of NUT in 293T lysates  
using NUT antibody.

