

P2Y6 Rabbit mAb


ABZ2111

Product Overview


Name	P2Y6 Rabbit mAb
Accession(Primary)	Q15077
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 36 kDa; Observed MW: 36 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; IP-1:20-1:50

Target information (Q15077)


Gene Name	P2RY6
About Gene Symbol	



Immunocytochemistry analysis of P2Y6 (green) in A549 using P2Y6 antibody, and DAPI (blue).



Western blot analysis of P2Y6 in K562, HeLa lysates using P2Y6 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using P2Y6 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.