


## PDGFR alpha (7H6) Mouse mAb

ABZ4080

### Product Overview

|                             |                             |
|-----------------------------|-----------------------------|
| <b>Name</b>                 | PDGFR alpha (7H6) Mouse mAb |
| <b>Accession(Primary)</b>   | P16234                      |
| <b>Host Species</b>         | Mouse                       |
| <b>Clonality</b>            | Monoclonal                  |
| <b>Isotype</b>              | IgG1                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat           |
| <b>Recommended Dilution</b> | IHC-1:50-1:100              |

### Target information (P16234)

**Gene Name** PDGFRA**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded rat Skin Tissue using PDGFR alpha (7H6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using PDGFR alpha (7H6) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.