

PDK2 Rabbit mAb

ABZ2107


Product Overview

Name	PDK2 Rabbit mAb
Accession(Primary)	Q15119
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 46 kDa; Observed MW: 46 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:200-1:1000; IP-1:20

Target information (Q15119)

Gene Name PDK2

About Gene Symbol



Western blot analysis of PDK2 in rat Brain lysates using PDK2 antibody.

