

PELP1 Rabbit mAb


ABZ1996

Product Overview


Name	PELP1 Rabbit mAb
Accession(Primary)	Q8IZL8
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 120 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; FC-1:200-1:500

Target information (Q8IZL8)


Gene Name	PELP1
About Gene Symbol	PELP1 (Proline, Glutamate And Leucine Rich Protein 1) is a membrane receptor involved in receptor-mediated signaling. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mec



Western blot analysis of PELP1 in HeLa, A549, MCF-7, 293T lysates using PELP1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PELP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of PELP1 (green)
in hela using PELP1 antibody, and DAPI(blue)