

PI3 Kinase p110 delta Rabbit pAb


ABZ1572

Product Overview


Name	PI3 Kinase p110 delta Rabbit pAb
Accession(Primary)	O00329
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 119 kDa; Observed MW: 119 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O00329)

Gene Name	PIK3CD
About Gene Symbol	



Western blot analysis of PI3 Kinase p110 delta in HepG2 lysates using PI3 Kinase p110 delta antibody.



Immunohistochemistry analysis of paraffin-embedded mouse brain using PI3 Kinase p110 delta antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.