

PKM2 (1F2) Mouse mAb


ABZ4167

Product Overview


Name	PKM2 (1F2) Mouse mAb
Accession(Primary)	P14618
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Species Reactivity	Human, Mouse, Cat
Calculated MW	Calculated MW: 58 kDa; Observed MW: 58 kDa
Recommended Dilution	WB-1:2000-1:8000; IF-1:50-1:200; IP: 1:20-1:50

Target information (P14618)


Gene Name	PKM
About Gene Symbol	PKM (Pyruvate Kinase M1/2) is a receptor tyrosine kinase involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of PKM2 in HCT116, COS7, K562, A549 and HeLa lysates using PKM2 antibody.



Immunofluorescence analysis of PKM2 (1F2) in 3T3 using PKM2 antibody.



Immunoprecipitation analysis of PKM2 (1F2) in HeLa lysates using PKM2 antibody.