

PP1C alpha (1C11) Mouse mAb

ABZ3930


Product Overview

Name	PP1C alpha (1C11) Mouse mAb
Accession(Primary)	P62136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IHC-1:50-1:100; IF-1:50-1:200


Target information (P62136)

Gene Name PPP1CA


About Gene Symbol



Western blot analysis of PPP1A in A431, Jurkat and HeLa whole lysates using PPP1A antibody.



Immunocytochemistry analysis of PP1C alpha (1C11) in Hela using PPP1A antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PPP1A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.