

PP1C alpha (1H1) Mouse mAb

ABZ3931


Product Overview

Name	PP1C alpha (1H1) Mouse mAb
Accession(Primary)	P62136
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:1000; IF-1:50-1:200


Target information (P62136)

Gene Name PPP1CA

About Gene Symbol



Western blot analysis of PPP1A in HepG2, MCF-7 and HeLa lysates using PPP1A antibody.



Immunocytochemistry analysis of PP1C alpha in HeLa cells using PP1C alpha antibody.