

Phospho-CDK1/2/3 (Thr14) Rabbit mAb


ABZ4532

Product Overview

Name	Phospho-CDK1/2/3 (Thr14) Rabbit mAb
Accession(Primary)	P06493
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 34 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:200; ICC/IF-1:100-1:500; IP-1:20-1:50

Target information (P06493)

Gene Name	CDK1/CDK2/CDK3
About Gene Symbol	



Western blot analysis of Phospho-
CDK1/CDK2/CDK3 (Thr14) in HeLa, CHO-K1 lysates
using Phospho-CDK1/2/3 (Thr14) antibody.