

Phospho-Cardiac Troponin I (Ser22/Ser23) Rabbit pAb


ABZ4480

Product Overview


Name	Phospho-Cardiac Troponin I (Ser22/Ser23) Rabbit pAb
Accession(Primary)	P19429
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 28 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P19429)


Gene Name	Tnni3
About Gene Symbol	




Western blot analysis of Phospho-Cardiac Troponin I (Ser22/Ser23) in SH-SY5Y lysates using Phospho-Cardiac Troponin I (Ser22/Ser23) antibody.



EnzymeLinked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phospho-peptide (Phospho-left) and NonPhospho-peptide (Phospho-right), using TNNI3 (Phospho-Ser22+Ser23) antibody.



Immunohistochemistry analysis of paraffin-embedded Human skeletal muscle using TNNI3 (Phospho-Ser22+Ser23) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Western blot analysis of Phospho-Cardiac Troponin I (Ser22/Ser23) in mouse heart lysates using Phospho-Cardiac Troponin I (Ser22/Ser23) antibody. The lane on the right is blocked with the Phospho-peptide.