

## Phospho-DDIT3 (Ser30) Rabbit pAb


ABZ4502

### Product Overview

|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-DDIT3 (Ser30) Rabbit pAb                                    |
| <b>Accession(Primary)</b>   | P35638                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse                                                        |
| <b>Calculated MW</b>        | Calculated MW: 19 kDa; Observed MW: 19 kDa                          |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:20000 |

### Target information (P35638)

|                          |                                                                                                                                                                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | DDIT3                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | DDIT3 (DNA Damage Inducible Transcript 3) is an epigenetic and nuclear regulatory protein involved in DNA damage response and repair. It is an important research target in epigenetic regulation, chromatin biology, and cancer research and is widely studied |



Western blot analysis of Phospho-CHOP (Ser30) in Jurkat lysates using Phospho-CHOP (Ser30) antibody.