

## Phospho-FAK (Tyr397) Rabbit mAb


ABZ4373

### Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Phospho-FAK (Tyr397) Rabbit mAb              |
| <b>Accession(Primary)</b>   | Q05397                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                            |
| <b>Calculated MW</b>        | Calculated MW: 119 kDa; Observed MW: 119 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:50-1:100           |

### Target information (Q05397)

|                          |                                                                                                                                                                                                                                                                                                |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | PTK2                                                                                                                                                                                                                                                                                           |
| <b>About Gene Symbol</b> | PTK2 (Protein Tyrosine Kinase 2) is a receptor tyrosine kinase involved in cell migration and invasion. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Phospho-FAK (Tyr397) in C6 lysates using Phospho-FAK (Tyr397) antibody.