

Phospho-FOXO1/3 (Ser322/S325) Rabbit pAb


ABZ4370

Product Overview

Name	Phospho-FOXO1/3 (Ser322/S325) Rabbit pAb
Accession(Primary)	Q12778
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 70 kDa; Observed MW: 97 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (Q12778)

Gene Name	FOXO1/FOXO3
About Gene Symbol	FOXO1 (Forkhead Box O1) is a transcription factor involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-FOXO1/3 (Ser322/S325) in 3T3 lysates using Phospho-FOXO1/3 (Ser322/S325) antibody.