

## Phospho-HDAC6 (Ser22) Rabbit pAb


ABZ4344

### Product Overview

|                             |                                                                     |
|-----------------------------|---------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-HDAC6 (Ser22) Rabbit pAb                                    |
| <b>Accession(Primary)</b>   | Q9UBN7                                                              |
| <b>Host Species</b>         | Rabbit                                                              |
| <b>Clonality</b>            | Polyclonal                                                          |
| <b>Isotype</b>              | IgG                                                                 |
| <b>Species Reactivity</b>   | Human, Mouse                                                        |
| <b>Calculated MW</b>        | Calculated MW: 131 kDa; Observed MW: 131 kDa                        |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:200-1:1000; ELISA-1:10000 |

### Target information (Q9UBN7)

|                          |                                                                                                                                                                                                                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | HDAC6                                                                                                                                                                                                                                                                                            |
| <b>About Gene Symbol</b> | HDAC6 (Histone Deacetylase 6) is a histone deacetylase involved in chromatin remodeling and epigenetics. It is an important research target in metabolic disease, cancer metabolism, and drug discovery and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of Phospho-HDAC6 (Ser22) in 3T3 lysates using Phospho-HDAC6 (Ser22) antibody.