

# Phospho-Histone H4 (Ser1) Rabbit pAb

ABZ4486


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Phospho-Histone H4 (Ser1) Rabbit pAb       |
| <b>Accession(Primary)</b>   | P62805                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 11 kDa; Observed MW: 17 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:2000                           |

## Target information (P62805)

**Gene Name** H4C1

**About Gene Symbol**



Western blot analysis of Phospho-Histone H4 (Ser1) in NIH/3T3 lysates using Phospho-STK4T387 antibody. NIH/3T3 were treated by CIP at 37°C for 1 hour or Nocodazole at 37°C for 20 hours.

