

# Phospho-IKB alpha (Ser32) Rabbit mAb


ABZ4437

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Phospho-IKB alpha (Ser32) Rabbit mAb       |
| <b>Accession(Primary)</b>   | P25963                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                          |
| <b>Calculated MW</b>        | Calculated MW: 36 kDa; Observed MW: 36 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IP-1:20                   |

## Target information (P25963)

|                          |                                                                                                                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | NFKBIA                                                                                                                                                                                                                                                                    |
| <b>About Gene Symbol</b> | NFKBIA (I-kappa-B-alpha) is a transcription factor involved in NF- $\kappa$ B signaling. It is an important research target in gene regulation, development, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets. |



Western blot analysis of IKB alpha (Phospho-Ser32) in 3T3, HeLa lysates using Phospho-IKB alpha (Ser32) antibody.