

Phospho-IKK alpha (Thr23) Rabbit pAb


ABZ4464

Product Overview


Name	Phospho-IKK alpha (Thr23) Rabbit pAb
Accession(Primary)	O15111
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 85 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:5000

Target information (O15111)

Gene Name	CHUK
About Gene Symbol	CHUK (Component Of Inhibitor Of Nuclear Factor Kappa B Kinase Complex) is a protein kinase and signaling regulator involved in NF- κ B signaling. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely st



Western blot analysis of Phospho-IKK alpha (Thr23) in 293 lysates using Phospho-IKK alpha (Thr23) antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma tissue using IKK α (Phospho-Thr23) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.