

# Phospho-KCNQ2/3/4/5 (Thr217/Thr246/Thr223/Thr251) Rabbit pAb


ABZ4485

## Product Overview


|                             |                                                                   |
|-----------------------------|-------------------------------------------------------------------|
| <b>Name</b>                 | Phospho-KCNQ2/3/4/5<br>(Thr217/Thr246/Thr223/Thr251) Rabbit pAb   |
| <b>Accession(Primary)</b>   | O43526                                                            |
| <b>Host Species</b>         | Rabbit                                                            |
| <b>Clonality</b>            | Polyclonal                                                        |
| <b>Isotype</b>              | IgG                                                               |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                                                 |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; IHC-1:100-1:300;<br>IF-1:50-1:200; ELISA-1:20000 |

## Target information (O43526)

|                          |                         |
|--------------------------|-------------------------|
| <b>Gene Name</b>         | KCNQ2/KCNQ3/KCNQ4/KCNQ5 |
| <b>About Gene Symbol</b> |                         |



EnzymeLinked Immunosorbent Assay (Phospho-ELISA) for Immunogen Phospho-peptide (Phospho-left) and NonPhospho-peptide (Phospho-right), using Kv7.3/KCNQ3 (Phospho-Thr24antibody



Immunohistochemical analysis of paraffin-embedded Human tonsils using Phospho-KCNQ2/3/4/5 (Thr217/Thr246/Thr223/Thr251) antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.