

# Phospho-MCM2 (Ser108) Rabbit mAb

ABZ4395


## Product Overview

|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Phospho-MCM2 (Ser108) Rabbit mAb             |
| <b>Accession(Primary)</b>   | P49736                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                            |
| <b>Calculated MW</b>        | Calculated MW: 102 kDa; Observed MW: 125 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; ICC/IF-1:20-1:50; IP-1:20   |


## Target information (P49736)

**Gene Name** MCM2


**About Gene Symbol**



Western blot analysis of Phospho-MCM2 (Ser108) in rat spleen lysates using Phospho-MCM2 (Ser108) antibody.



Western blot analysis of Phospho-MCM2 (Ser108) in K562, 3T3, HeLa lysates using Phospho-MCM2 (Ser108) antibody.



Immunocytochemistry analysis of MCM2 (Phospho-Ser108)(green) in HeLa using MCM2 (Phospho-Ser108) antibody, and DAPI(blue)