

Phospho-MCM2 (Ser27) Rabbit mAb


ABZ4396

Product Overview


Name	Phospho-MCM2 (Ser27) Rabbit mAb
Accession(Primary)	P49736
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 102 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:100-1:200; IP-1:10-1:100

Target information (P49736)


Gene Name	MCM2
About Gene Symbol	



Immunocytochemistry analysis of Phospho-MCM2 (Ser27) (green) in HeLa using Phospho-MCM2 (Ser27) antibody, and DAPI (blue).



Western blot analysis of MCM2 (Phospho-Ser27) in HeLa, A549, HL-60, U2OS lysates using MCM2 (Phospho-Ser27) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using MCM2 (Phospho-Ser27) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.