

Phospho-PDPK1 (Ser241) Rabbit pAb

ABZ4583

Product Overview

Name	Phospho-PDPK1 (Ser241) Rabbit pAb
Accession(Primary)	O15530
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 63 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (O15530)

Gene Name	PDPK1
About Gene Symbol	PDPK1 (3-Phosphoinositide Dependent Protein Kinase 1) is a protein kinase and signaling regulator involved in cell signaling and proliferation. It is an important research target in cancer, cell signaling, and therapeutic target research and is widely stu

