

## Phospho-PML Protein (Ser518) Rabbit mAb


ABZ4571

### Product Overview


|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | Phospho-PML Protein (Ser518) Rabbit mAb        |
| <b>Accession(Primary)</b>   | P29590                                         |
| <b>Host Species</b>         | Rabbit                                         |
| <b>Clonality</b>            | Monoclonal                                     |
| <b>Isotype</b>              | IgG                                            |
| <b>Species Reactivity</b>   | Human, Mouse                                   |
| <b>Calculated MW</b>        | Calculated MW: 98 kDa; Observed MW: 98,117 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:100-1:200            |

### Target information (P29590)

|                          |     |
|--------------------------|-----|
| <b>Gene Name</b>         | PML |
| <b>About Gene Symbol</b> |     |



Western blot analysis of Phospho-PML Protein (Ser518) in A549, HL-60, U2OS lysates using Phospho-PML Protein (Ser518) antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PML Protein (Phospho-Ser518) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.