

Phospho-Phospholamban (Ser16) Rabbit mAb

ABZ4436

Product Overview

Name	Phospho-Phospholamban (Ser16) Rabbit mAb
Accession(Primary)	P26678
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 6 kDa; Observed MW: 12,24 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (P26678)

Gene Name

PLN

About Gene Symbol


