

Phospho-RSK4 (Ser232) Rabbit mAb


ABZ4511

Product Overview



Name	Phospho-RSK4 (Ser232) Rabbit mAb
Accession(Primary)	Q9UK32
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:1000-1:5000; IHC-P-1:50-1:100; ICC/IF-1:50-1:100

Target information (Q9UK32)

Gene Name	RPS6KA6
About Gene Symbol	



Immunocytochemistry analysis of Phospho-RSK4 (Ser232) (green) in LNCap using Phospho-RSK4 (Ser232) antibody, and DAPI (blue).



Western blot analysis of Phospho-RSK4 (Ser232) in HeLa lysates using Phospho-RSK4 (Ser232) antibody.