

## Phospho-Rb (Thr356) Rabbit mAb


ABZ4587

### Product Overview


|                             |                                              |
|-----------------------------|----------------------------------------------|
| <b>Name</b>                 | Phospho-Rb (Thr356) Rabbit mAb               |
| <b>Accession(Primary)</b>   | P06400                                       |
| <b>Host Species</b>         | Rabbit                                       |
| <b>Clonality</b>            | Monoclonal                                   |
| <b>Isotype</b>              | IgG                                          |
| <b>Species Reactivity</b>   | Human                                        |
| <b>Calculated MW</b>        | Calculated MW: 106 kDa; Observed MW: 110 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000; IHC-P-1:100-1:200; IP-1:20  |

### Target information (P06400)

|                          |                                                                                                                                                                                                                                                                      |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>         | RB1                                                                                                                                                                                                                                                                  |
| <b>About Gene Symbol</b> | <p>RB1 (RB Transcriptional Corepressor 1) is a protein phosphatase involved in gene transcription and regulation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.</p> |



Western blot analysis of Rb (Phospho-Thr356) in HL-60, U2OS, C6 lysates using Phospho-Rb (Thr356) antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using Rb (Phospho-Thr356) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.