

## Phospho-TrkB (Tyr705) Rabbit mAb


ABZ4351

### Product Overview

|                             |                                                |
|-----------------------------|------------------------------------------------|
| <b>Name</b>                 | Phospho-TrkB (Tyr705) Rabbit mAb               |
| <b>Accession(Primary)</b>   | Q16620                                         |
| <b>Host Species</b>         | Rabbit                                         |
| <b>Clonality</b>            | Monoclonal                                     |
| <b>Isotype</b>              | IgG                                            |
| <b>Species Reactivity</b>   | Human, Mouse, Rat                              |
| <b>Calculated MW</b>        | Calculated MW: 92 kDa; Observed MW: 140,90 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:1000                                |

### Target information (Q16620)

|                          |       |
|--------------------------|-------|
| <b>Gene Name</b>         | NTRK2 |
| <b>About Gene Symbol</b> |       |



Western blot analysis of Phospho-TrkB (Tyr705) in HeLa, A549, C6 lysates using Phospho-TrkB (Tyr705) antibody.