

Phospho-delta 1 Catenin (Thr916) Rabbit mAb


ABZ4558

Product Overview


Name	Phospho-delta 1 Catenin (Thr916) Rabbit mAb
Accession(Primary)	O60716
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 108 kDa; Observed MW: 100 kDa
Recommended Dilution	WB-1:2000-1:5000; IHC-P-1:100-1:200; ICC/IF-1:50-1:100

Target information (O60716)


Gene Name	CTNND1
About Gene Symbol	CTNND1 (Catenin Delta 1) is a cell junction/adhesion protein involved in cytoskeletal organization. It is an important research target in cell biology, cancer metastasis, and cytoskeletal research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.



Western blot analysis of Phospho-delta 1 Catenin/CAS (Thr916) in MCF-7, U251, C2C12 lysates using Phospho-delta 1 Catenin (Thr916) antibody.



Immunocytochemistry analysis of delta 1 Catenin/CAS (Phospho-Thr916) (green) in HeLa using delta 1 Catenin/CAS (Phospho-Thr916) antibody, and DAPI (blue)



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Phospho-delta 1 Catenin/CAS (Thr916) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.