


# Prostate Secretory Protein Rabbit mAb

ABZ5190

## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Prostate Secretory Protein Rabbit mAb      |
| <b>Accession(Primary)</b>   | P08118                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Monoclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 13 kDa; Observed MW: 13 kDa |
| <b>Recommended Dilution</b> | WB-1:1000-1:5000; IHC-P-1:200-1:1000       |

## Target information (P08118)

**Gene Name** MSMB**About Gene Symbol**

Western blot analysis of Prostate Secretory Protein/PSP in Human prostate cancer lysates using Prostate Secretory Protein antibody.

