


RAB45 Rabbit mAb


ABZ1997

Product Overview


Name	RAB45 Rabbit mAb
Accession(Primary)	Q8IZ41
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 83 kDa; Observed MW: 83 kDa
Recommended Dilution	WB-1:1000-1:5000; ICC/IF-1:20-1:50; FC-1:20-1:50

Target information (Q8IZ41)

Gene Name RASEF**About Gene Symbol**



Western blot analysis of RAB45 in Jurkat lysates using RAB45 antibody.



Immunocytochemistry analysis of RAB45 (green) in HeLa using RAB45 antibody, and DAPI (blue)