


RAGE Rabbit mAb

ABZ2109

Product Overview

Name	RAGE Rabbit mAb
Accession(Primary)	Q15109
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Species Reactivity	Mouse, Rat
Calculated MW	Calculated MW: 43 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:1000-1:5000

Target information (Q15109)

Gene Name AGER**About Gene Symbol**

Western blot analysis of RAGE in mouse lung
lysates using RAGE antibody.

