

AHSP Human

Alpha Hemoglobin Stabilizing Protein Human Recombinant

RCP0035

Product Overview

Name	AHSP Human
Accession(Primary)	Q9NZD4
Immunogen	Anti-human AICDA mAb, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human AICDA amino acids 1-54 purified from E. coli.
Description	Alpha Hemoglobin Stabilizing Protein Human Recombinant

Target information (Q9NZD4)

Gene Name

About Gene Symbol