

ASF1A Human

ASF1 Anti-Silencing Function 1 Homolog A Human Recombinant

RCP0126

Product Overview

Name	ASF1A Human
Accession(Primary)	Q9Y294
Immunogen	Anti-human ASF1B mAb, clone PAT1D4AT, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ASF1B protein (1-202 amino acids) purified from E. coli.
Description	ASF1 Anti-Silencing Function 1 Homolog A Human Recombinant

Target information (Q9Y294)

Gene Name

About Gene Symbol