

BAG1 Human

BCL2-Associated Athanogene 1 Human Recombinant

RCP0165

Product Overview

Name	BAG1 Human
Accession(Primary)	Q99933
Immunogen	Anti-human BAG2 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human BAG2 protein 1-211 amino acids purified from E. coli.
Description	BCL2-Associated Athanogene 1 Human Recombinant

Target information (Q99933)

Gene Name

About Gene Symbol

BAG1 (BAG Cochaperone 1) is a molecular chaperone involved in apoptosis. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.