

## BAIAP2 Human

*BAI1-Associated Protein 2 Human Recombinant*

RCP0168

### Product Overview

|                           |                                                                                                                                                                                                 |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>               | BAIAP2 Human                                                                                                                                                                                    |
| <b>Accession(Primary)</b> | Q9UQB8                                                                                                                                                                                          |
| <b>Immunogen</b>          | Anti-human BAK1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human BAK1 amino acids 29-187 purified from E. coli. |
| <b>Description</b>        | BAI1-Associated Protein 2 Human Recombinant                                                                                                                                                     |

### Target information (Q9UQB8)

**Gene Name**

**About Gene Symbol**