

BCAR1 Human

Breast Cancer Anti-Estrogen Resistance 1 Human Recombinant

RCP0178

Product Overview

Name	BCAR1 Human
Accession(Primary)	P56945
Immunogen	Anti-human BCAT1 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with recombinant human BCAT1 amino acids 1-386 purified from E. coli.
Description	Breast Cancer Anti-Estrogen Resistance 1 Human Recombinant

Target information (P56945)

Gene Name	
About Gene Symbol	BCAR1 (BCAR1 Scaffold Protein, Cas Family Member) is a signaling scaffold protein involved in cell signaling and proliferation. It is an important research target in cancer research and is widely studied to elucidate cellular mechanisms and identify therapeutic targets.