

BTN3A1 Human

Butyrophilin Subfamily 3 Member A1 Human Recombinant

RCP0242

Product Overview

Name	BTN3A1 Human
Accession(Primary)	O00481
Immunogen	Anti-human BUB3 mAb, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human BUB3 protein 1-328 amino acids purified from E. coli.
Description	Butyrophilin Subfamily 3 Member A1 Human Recombinant

Target information (O00481)

Gene Name

About Gene Symbol